


# HPLC PAL Systems

## Front End Automation Systems for Liquid Chromatography



### Drug Discovery LCMS

**PAL** Syringe only concept, no tubing in sample path for transparent sample injection

**PAL** Reliability and ruggedness for unattended 24 hour/day chromatography

**PAL** High sample capacity combined with fast injection cycles

**PAL** Temperature controlled sample storage from 4iC up to 70iC

**PAL** Third party instrument drivers for Micromass Masslynx<sup>a</sup> PE SCIEX Analyst<sup>a</sup> Agilent ChemStation<sup>a</sup> TMQ Xcalibur<sup>a</sup>

High Throughput Screening

Combinatorial Chemistry

Flow Injection Analysis


Preclinical Research

Pharmaco-kinetic studies

# HTS PAL

## Reliability is standard

Automated pretreatment operations combined with precise injection of sample sequences is certainly a key factor in increased laboratory productivity. Today, requirements in a sample injection system include high precision and accuracy, random access to individual samples and task-oriented programming along with options and accessories for additional sample handling routines. CTC Analytics' experience in the design and manufacture of automated GC and LC injection systems is reflected in their range of the PAL series of Prep and Load Systems. In all PAL instruments the latest electronic and mechanical developments are combined with a unique, stepper motor driven XYZ design. This combination, together with complete protocols for IQ/OQ instrument qualification, ensures you precise, fast sample handling for dependable analysis results.



## Specifications HTS PAL

### Syringe size:

100µl standard  
(Optional syringe sizes 10, 25, 250, 1000, 2500, 5000µl)

### Injection volume range:

0.1µl - 5000µl

### Injection speed:

0.01µl/sec. - 250µl/sec

### LC Injector:

1 electrically actuated fast switching 6-port valve  
(optional second injection valve available)

### Sample capacity:

up to 600  
 2ml standard vials  
 1 96 0ml or 20ml vials  
 deepwell plates ( 96 or 384 wells)  
 standard microplates ( 96 or 384 wells)

### Needle, Syringe and Valve cleaning:

Fast Wash Station for 2 different solvents

### Remote control:

Cycle Composer Windows 9x/NT4/2000 based PC software

### Electrical control:

2 RS 232C ports / 3 TTL Input / 3 TTL Output /  
2 Opto Coupler Input / 2 Relay Output

### Dimensions:

L: 828mm D:385mm H: 648mm

### Weight:

10kg (without accessories)

### HTS PAL Options:

Thermostatted Trayholders (4iC ÷ 70iC)  
Thermostatted Microplate Stacks (4iC ÷ 40iC)  
4-port and 10-port electrically actuated fast switching injection valves  
Barcode Reader  
PAL Dilutor for Dilutions / Derivatisations



Temperature controlled sample storage (4iC-40iC) for Micro- or Deepwell plates



Additional capabilities using the PAL Dilutor



Barcode reading ensures positive sample identification

# HTC PAL

## Advanced performance

The PAL series of LC sample loaders provide outstanding performance and maximum flexibility for your HPLC sample processing system. Sample capacity of up to 24 microplates within 50cm benchspace are unmatched in the industry. Temperature controlled sample storage makes it easy to cool down samples to prevent degradation or heat samples for kinetic studies. Using the optional 4-port or 10-port injection valves enables complex column switching procedures such as pre-column sample enrichment or internal loop configurations.

## Maintenance

CTCO's PAL instrument design provides worry-free operation and low maintenance costs. An open architecture for easy access to the injection valves, sample trays and syringes, guarantees the quick exchange of sample loops, trays and different syringe sizes. All PAL liquid handling systems are equipped with Flash EEPROM technology, for problem-free update of every installed unit with the newest enhancements and capabilities.

## Specifications HTC PAL

### Syringe size:

100µl standard  
(Optional syringe sizes 10, 25, 250, 1000, 2500, 5000µl)

### Injection volume range:

0.1µl - 5000µl

### Injection speed:

0.01µl/sec. - 250µl/sec

### LC Injector:

1 electrically actuated fast switching 6-port valve

### Sample capacity:

up to 400ml micro vials

2ml vials

1 64 0ml or 20ml vials

deepwell ( 96 or 384 wells)

standard microplates ( 96 or 384 wells)

### Needle, Syringe and Valve cleaning

Fast Wash Station for 2 different solvents

### Remote control

Cycle Composer Windows 9x/NT4/2000 based PC software

### Electrical control

2 RS 232C ports / 3 TTL Input /

1 Opto Coupler Input / 2 Relay Output

### Dimensions

L: 534mm D:385mm H: 648mm

### Weight

8kg (without accessories)

### HTC PAL Options


Thermostatted Trayholders (4iC ÷ 70iC)

Thermostatted Microplate Stacks (4iC ÷ 40iC)

4-port and 10-port electrically actuated fast switching injection valves

Barcode Reader


Diluent Reservoirs for Dilutions / Derivatisations



The HTC/HTS PAL automatically opens a stack drawer and aspirates the sample from the microplate



and injects it directly into the fast switching LC valve



Fast Wash Station for high throughput applications



**Syringe only technology**

The syringe only concept of the PAL series of LC sample loaders combines the manual sample injection procedure with the precision and throughput of a robotic liquid handling system. Samples are aspirated with a conventional liquid syringe and are transferred directly into the sample loop. No error prone teflon tubing or transfer lines are involved during the sample loading process. Before and after each injection the syringe and syringe needle (needle inside and outside) are rinsed in a dedicated Fast Wash Station containing two different cleaning solvents. This unique design virtually eliminates sample carryover and adsorption problems.

**Large sample capacity combined with fast cycle time**


You expect an autosampler to provide not only precision, but also reliable, fast injections to enhance laboratory productivity. CTC's unique XYZ mechanism allows for fast sampling (2-3 injections per minute) or rapid vial to vial transfers. Fast injection cycles coupled with a large sample capacity, makes the PAL series of LC sample loaders ideal for high throughput screening, QA/QC applications or LCMS Flow Injection Analysis (FIA).



Choose from a wide variety of syringe sizes

**Specifications LC PAL**

<b>Syringe size:</b> 100µl standard (Optional syringe sizes 10, 25, 250, 1000, 2500, 5000µl)
<b>Injection volume range:</b> 1.0µl - 5000µl
<b>Injection speed:</b> 0.01µl/sec. - 250µl/sec
<b>LC Injector</b> 1 electrically actuated fast switching 6-port valve
<b>Sample capacity</b> up to 216 2ml vials 4 deepwell or standard microplates ( 96 or 384 wells)
<b>Needle, Syringe and Valve cleaning</b> Wash Station for 2 different solvents
<b>Remote control</b> Cycle Composer Windows 9x/NT4/2000 based PC software
<b>Electrical control</b> 2 RS 232C ports / 3 TTL Input / 1 Opto Coupler Input / 2 Relay Output
<b>Dimensions</b> L: 534mm D:385mm H: 545mm
<b>Weight</b> 8kg (without accessories)



Log book with protocols for IQ/OQ instrument qualification



4- 6- and 10 port Injection Valves used for different applications



# Cycle Composer

## Intelligent Automation

The Windows 9x/NT4/2000 software Cycle Composer provides remote control for the PAL series of LC sample loaders. The Cycle Composer software allows the operator to easily setup, edit and run HTS / HTC and LC PAL methods for even very complex "Prep and Load" applications.

## Easy to use

The Cycle Composer has a graphical workspace editor that allows the user easily to configure the workspace, accepting various microplates (96 and 384 wells), 1ml microvials or 2ml standard vials. These various formats may be matched or mixed within the workspace. Samples can be selected serially or randomly for processing. The "Point and click" operation quickly directs the user through programming steps that configure instrument setup, methods and sample lists. For routine daily use, no special programming skills are needed. The Cycle Composer contains a library of common liquid handling procedures including sample transfer, reagent/standard addition, mixing and dilution steps.

## Customize your HPLC PAL

As with the PAL hardware concept the Cycle Composer is already prepared for individual application requirements. Additional flexibility can be assigned to sample prep procedures by using the powerful Cycle Composer macro language, which expands the application range of every HTS / HTC or LC PAL System.

## Third party data acquisition drivers:

For "single keyboard" operation, various third party data acquisition drivers are available. This enables the complete integration of all Cycle Composer capabilities into one single LC processing system.




## Cycle Composer specifications:

### Min. PC requirements:


- Pentium processor >90MHz
- 32 MB RAM
- 1 serial interface
- 1 3.5" floppy drive
- 1 CD ROM drive
- installed Windows 9x/NT 4.0/2000 or higher
- 10 MB free Harddisk space

### Compatible PAL Instruments:


HTS PAL, HTC PAL, LC PAL  
(including PAL Firmware 2.2 or higher)



Method editor: Setup and edit customized injection cycles



Sample list: Setup and run injection sequences



Trays: Graphical sample list generator



## General Specifications

### System type:

XYZ robot with syringe only concept, no tubing in sample path

### Local user interface:

Control panel with 4 function keys, graphical LCD display, unique scroll knob for teach functions

### Remote control:

Cycle Composer control software Windows 9x/NT4/2000

### Maintenance:

Accessibility to all maintenance parts from front  
Preventative maintenance kits available

### Power Requirements:

97.5 ± 264 VAC, 45-66 Hz

### Environment:

10°C - 40°C constant temperature,  
< 80% humidity (non condensing)

### Injection Parameter Control

- Injection volume
- Pre - injection syringe wash strokes for two different solvents
- Post - injection syringe wash strokes for two different solvents
- Pre injection syringe wash strokes with sample
- Plunger speed used to aspirate / eject sample
- Air gap after sample aspiration
- Number of filling strokes to aspirate sample
- Delay time between sample pull-up and ejection
- Injector selection used for injection
- Plunger speed used during sample injection
- Delay time prior and after sample injection
- Valve cleaning wash strokes

Specifications are subject to change without notice

## ...other PAL sample injection systems

COMBI PAL Front End Automation  
for Gaschromatography



Distributed by:

CTC Analytics has dedicated the last 10 years to the continued development and high reliability of advanced sample injection technology. To learn more about the unique PAL Series of LC/LCMS sample handling systems or any of our GC/GCMS sample injectors contact your CTC Analytics distributor.

CTC Analytics AG acknowledges all tradenames and trademarks used as the property of their respective owners.

CTC Analytics AG  
Industriestrasse 20  
CH-4222 Zwingen  
Switzerland  
Tel: +41 61 7658100  
Fax: +41 61 7658199  
E-mail: info@ctc.ch  
Web: www.ctc.ch



...CTC ANALYTICS

Where design meets performance